

# MANUAL

## **BIX-LV19 Sterilized Training Manikin**



**Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.**

**Ada Med Supply Limited**

## **BIX-LV19 Sterilized Training Manikin**

### **一、 Features Of Model:**

The model is adult male upper torso, with the anatomical structure of the chest, realistic shape, clearly marked characteristics. This model is made of imported PVC plastic materials, use stainless steel molds, through injection pressure of high temperatures to take its shape.

### **二、 Function:**

1. With good flexibility, elasticity and lifelike texture, the model provides a good training aid for aseptic skill.
2. Human body surface accurately signs: thoracic, clavicle, chest Angle rib bones, ribs, clearance, ventral horn, xiphoid process, rib on the edge of the bow, before the ilium, the iliac spine and pubic bone, etc.
3. Can do surgery area disinfection and spread towels operation training

### **三、 Maintenance And Warranty:**

1. Keep the model clean, after the operation and training end, the model should be wiped clean.
2. When long time do not use the model, please wipe clean the model, packing it and stored in a ventilated, cool and dry place to prolong life.

**Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.**

**Ada Med Supply Limited**

---

Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd. Building 3,  
No.1288 kungang Road, Xiaokunshan Town, Songjiang District, Shanghai, China

TEL: +86-13383897707

FAX: 021—57722933

Post code: 201612

E-mail: [edith@adahealthy.com](mailto:edith@adahealthy.com)

Website: [www.adahealthy.com](http://www.adahealthy.com)