

# MANUAL

## **BIX/F132 Neonatal umbilical cord care model**



**Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.**

**Ada Med Supply Limited**

## BIX/F132 Neonatal umbilical cord care model

Neonatal umbilical cord care model is made of imported PVC material, which is characterized by lifelike image, real operation, reasonable structure and durability. It can be used as baby bath, changing clothes and diapers and umbilical care, eye care and ear care. Guidance on nursing training and teaching.

### **Features:**

1. Full body Neonatal model, head, limbs are flexible.
2. Can do Neonatal breastfeeding, bathing, dressing, changing diapers and other basic care operations. Can do neonatal umbilical cord ligation and nursing operation training. Real feeling, umbilical cord with a certain length can be used for many times.
3. Realistic umbilical, umbilical vein and umbilical artery inside maintenance.

### **Maintenance :**

1. After the end of the operation training, especially when not in use for a long time, the model should be cleaned clean, safe and stored in a cool dry place.
2. If the model surface is contaminated, use a soft cleaning cloth dipped in detergent or soapy water to gently wipe away.
3. The model can not touch benzene, xylene, acetone, chloroform, gasoline and other organic solvents.

**Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.**

**Ada Med Supply Limited**

---

Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd. Building 3,  
No.1288 kungang Road, Xiaokunshan Town, Songjiang District, Shanghai, China

Tel: +86-13383897707

Fax: 021—57722933

Post code: 201612

E-mail: [edith@adahealthy.com](mailto:edith@adahealthy.com)

Website: [www.adahealthy.com](http://www.adahealthy.com)