

MANUAL

BIX/F7B Advanced breast diagnosis and palpation model



Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.

Ada Med Supply Limited

BIX/F7B Advanced breast diagnosis and palpation mode

I.Features

Typical symptoms of various common breast tumor, designed specifically for female breast clinical diagnosis and self-examination.

There are more than 10 kinds of different lesions will meet the needs of teaching.

The model provides the following pathological changes:

Nodules: Texture hard, surface is not smoothly, can be considered as malignant tumor ; texture relatively soft, surface smoothing, can be considered as benign tumor.

lymphatic metastasis : Hard lymph nodules can be hit on lymph nodules

The change of the nipple, inverted nipple; Nipple has burst and bloody fluid overflow.

The change of the skin: Skin sag; appearance looks like orange peel; Inflammatory breast cancer.

II.Application notice

- 1.Avoid to let the model long time exposure in the sunlight to prevent deformation and discoloration.
2. When the model is polluted, please use cotton gauze or cotton swab with diluted detergent to wash.
3. Try to avoid the model fall or collision, and any scratch of sharp objects on model or cutting and violent treatment will shorten service life.

Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.

Ada Med Supply Limited

Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd. Building 3,
No.1288 kungang Road, Xiaokunshan Town, Songjiang District, Shanghai, China

Tel: +86-13383897707

Fax: 021—57722933

Post code: 201612

E-mail: edith@adahealthy.com

Website: www.adahealthy.com