

# MANUAL

## **BIX/F14 Advanced Breast Examination Model**



**Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.**

**Ada Med Supply Limited**

## BIX/F14 Advanced Breast Examination Model

### Main functions of the product:

- A natural and realistic model shows the breasts of women. It can be used for training and practice of breast examination.

The pictures and text content of this product are provided by Shanghai Zhineng Medical Model Manufacturing Center and all rights of interpretation are assumed. If the reproduction or illegal publicity is carried out without consent, our company will pursue legal responsibility.

- The dress-style structure design can better carry out self-examination of breast training.

- Various common breast tumors are buried in various parts of the breast:

Hard texture and rough surface can be regarded as malignant tumors

Relatively soft texture and smooth surface can be regarded as benign tumors

Light texture hard palate lymph nodes

Strip-shaped lobular hyperplasia.

- Changes in nipples and skin:

Nipples and skin depressions

Orange peel appearance.

### Precautions

1. Avoid long-term exposure and sun exposure to the model to prevent deformation and discoloration.
2. When the model is contaminated, it can be wiped with gauze or cotton wool dipped in diluted detergent.
3. Be careful to avoid the model from falling or bumping. Any sharp objects scratching, cutting or violent treatment of the model will shorten the service life.

**Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.**

**Ada Med Supply Limited**

---

Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd. Building 3,

No.1288 kungang Road, Xiaokunshan Town, Songjiang District, Shanghai, China

Tel: +86-13383897707

Fax: 021—57722933

Post code: 201612

E-mail: [edith@adahealthy.com](mailto:edith@adahealthy.com)

Website: [www.adahealthy.com](http://www.adahealthy.com)