

# MANUAL

## BIX-LV11 ADVANCED TROCAR TRAINING MODEL



**Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.**

**Ada Med Supply Limited**

## **BIX-LV11 ADVANCED TROCAR TRAINING MODEL**

### **1.Features:**

The advanced trocar training model is a combination of commonly used clinical surgical procedures for laparoscopic surgery operation training, and is designed with reference to relevant teaching materials. It is the best practice way to improve surgical technique. Each model can be repeated for many technical exercises. The module can be replaced after it is discarded. It is one of the most commonly used teaching aids for clinical surgical operation training. Features such as lifelike shape, real operation, reasonable structure, firmness and durability, and non-deformation during disinfection and cleaning.

### **2.Functions:**

1. Transparent abdominal wall trocar insertion and puncture practice
2. Clear anatomical level, good elasticity of skin tissue, real hand feeling, realistic blocking feeling when the needle passes through the skin, and a sense of failure after the needle passes through the skin
3. Through the transparent seat, the whole process of trocar insertion can be observed
4. Various types of trocars can be used for insertion

### **3.Application:**

1. Clinical teaching demonstration and practical operation training of students in medical universities, nursing schools and medical colleges
2. Continuing education for hospital medical and nursing staff, clinical teaching practice operation training

General training of clinical medicine in primary health units

#### **4.Maintenance and warranty:**

1. After the operation and training of the model is over, it can be scrubbed and cleaned with soap and water. Do not expose the model to direct sunlight.
2. When not in use for a long time, the model should be wiped clean and stored in a cool and dry place to prolong its service life.

**Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.**

**Ada Med Supply Limited**

---

Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd. Building 3,  
No.1288 kungang Road, Xiaokunshan Town, Songjiang District, Shanghai, China

TEL: +86-13383897707

FAX: 021—57722933

Post code: 201612

E-mail: [edith@adahealthy.com](mailto:edith@adahealthy.com)

Website: [www.adahealthy.com](http://www.adahealthy.com)