

MANUAL

BIX-LV41 Advanced Ear Examination Model



Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.

Ada Med Supply Limited

BIX-LV41 Advanced Ear Examination Model

一、Advanced Ear Examination Model

This product is a multifunctional model for training ear examination and treatment.

The model accurately simulates the structure of the outer ear and middle ear, and can simulate ears with different diseases through simple disassembly and installation. The model is made in a life-size. By replacing the ear attachments with colored lesions, the trainers are provided with a variety of real scenes for inspection, training and assessment.

二、Main Features

- Standard ear examination position;
- Precise anatomical structure of auricle, external auditory canal and tympanic membrane;
- Ear lesions examination via an otoscope;

三、Ear Lesions Simulators Components:

1. Normal tympanic membrane;
2. Inverted tympanic membrane;
3. Small perforation of tympanic membrane;
4. Total perforation of tympanic membrane;
5. Traumatic tympanic membrane perforation;
6. Dry posterior central perforation;
7. Tympanic membrane incision and tube placement;
8. Bullous tympanitis;
9. Herpes blisters on the tympanic membrane;
10. Tympanic membrane sclerosis;
11. Tympanosclerosis Crescent sclerosis;

- 12. Exudation of serous otitis media;
- 13. Early congestion of acute otitis media;
- 14. Acute otitis media;
- 15. Suppurative otitis media;
- 16. Chronic suppurative otitis media;
- 17. Cholesteatoma;
- 18. Foreign body in the ear;
- 19. Cleaning in the ear;

Note: Cerumen cleaning operation is available and ear lesions simulators components are replaceable.

Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.

Ada Med Supply Limited

Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd. Building 3,
No.1288 kungang Road, Xiaokunshan Town, Songjiang District, Shanghai, China

TEL: +86-13383897707

FAX: 021—57722933

Post code: 201612

E-mail: edith@adahealthy.com

Website: www.adahealthy.com