

# MANUAL

## **BIX/H140 Senior infant nursing model**



**Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.**

**Ada Med Supply Limited**

## **BIX/H140 Senior infant nursing model**

### **Function characteristics:**

1. Baby head and neck and limbs can move freely, baby boy/girl chest skin can be replaced

Infant head L venipuncture

2. Vein injection, infusion (blood) piercing practice of baby right leg great saphenous vein, small saphenous

3. Baby umbilical cord vein puncture

4. Baby mouth nasal cannula nursing

5. Gastric lavage, enema, urethral catheterization

6. Colostomy drainage

7. Deltoid, buttock muscle injection

8. Overall care, bath, lactation, change clothes, change diapers

### **Application range**

1. Medical university, college of nursing, health school, the teaching of clinical medical colleges and universities teaching and student practical training.

2. Hospital medical and nursing staff to continue education in clinical teaching practice training.

3. Health units at basic level training of clinical popularization.

### **Maintenance warranty**

1. When model with besmirch, please clean it with neutral detergent, detergent, detergent, alcohol, etc.).

2. Avoid sharp objects cut on the model.

3. Avoid long-term sun exposure.

**Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.**

**Ada Med Supply Limited**

---

Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd. Building 3,  
No.1288 kungang Road, Xiaokunshan Town, Songjiang District, Shanghai, China

Tel: +86-13383897707

Fax: 021—57722933

Post code: 201612

E-mail: [edith@adahealthy.com](mailto:edith@adahealthy.com)

Website: [www.adahealthy.com](http://www.adahealthy.com)